

Technical Data Sheet

Product description: Flexible gingiva mask material for the C&B technique

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical term
1	polyvinylmethyilsiloxane
2	polymethylhydrogensiloxane
3	silicic acid
4	fillers
5	pigments
6	platinum-catalyst

Properties:

Color: pink (gingiva colored)

Processing time: 1 min 30 s $\pm$ 10 s at 23 °C

Deflasking time: 5 min $\pm$ 30 s at 37 °C

Final hardness: 60 $\pm$ 2 Shore A (DIN ISO 48-4)

Consistency: light-bodied - Type 3 (DIN EN ISO 4823)

Elastic recovery: > 99.5 % (DIN EN ISO 4823)

Linear dimensional change: < 0.5 % (DIN EN ISO 4823)