

# GumQuick implant



## Technical Data Sheet

Product description: Hard/flexible gingiva mask material for the implantology

Manufacturer: Dreve Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

### Composition:

| Pos | chemical term              |
|-----|----------------------------|
| 1   | polyvinylmethyilsiloxane   |
| 2   | polymethylhydrogensiloxane |
| 3   | silicic acid               |
| 4   | fillers                    |
| 5   | pigments                   |
| 6   | platinum-catalyst          |

### Properties:

Color: pink (gingiva colored)

Processing time: 1 min 30 s  $\pm$  10 s at 23 °C

Deflasking time: 5 min  $\pm$  30 s at 37 °C

Final hardness: 73  $\pm$  2 Shore A (DIN ISO 48-4)

Consistency: light-bodied - Type 3 (DIN EN ISO 4823)

Elastic recovery: > 99.5 % (DIN EN ISO 4823)

Linear dimensional change: < 0.5 % (DIN EN ISO 4823)

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@dreve.de

www.dreve.com