

## Angebot #442140

	Artikel	Stückpreis
 <p>RECIPROC Gutta-Percha ISO 050 • Precise fit to canals shaped with RECIPROC® NiTi instruments • All cold and warm obturation techniques • Phase with low melting point • Injection moulded</p> <p>Endo Easy Efficient®</p>	<p><b>Packung 60 Stück ISO 050</b></p> <p>Artikelnummer: 22201</p> <p><b>ISO Größe</b> <input type="text" value="050"/></p> <p><b>Verpackungseinheit</b> <input type="text" value="Packung"/></p> <p><b>Verpackungsinhalt</b> <input type="text" value="60 Stifte"/></p>	<p><b>20,21 €* </b></p>